
Configuration Instructions


Team 7 - CS 157A sec 03 - project 3

Requirements: Javascript, Node.js, mySQL v8.0, express, ejs


- A. Create database model
 - 1. Install and run mySQL v8.0
 - 2. Run mySQL Workbench



3. Select and connect to a connection as root or any user that can modify database



4. Open file "insert data.sql" and run all the queries




B. Set up to run program

- 1. Open Command Prompt on Windows or any other console.
- 2. Run npm install express ejs mysql.
- 3. Locate to where you save the zip file and unzip that file.
- 4. Run node app.js on command to start the webpage

```
C:\Users\ngoth\Documents\git hub\team-7-topic-3-cs157a-3\DifferentAproach>nodemon app.js
[nodemon] 2.0.1
[nodemon] to restart at any time, enter `rs`
[nodemon] watching dir(s): *.*
[nodemon] watching extensions: js,mjs,json
[nodemon] starting `node app.js`
Server running on port: 5000
```

- C. Interview Webpage Instruction
 - 1. Connection page open any browsers and access to
 - ① localhost:5000
 - 2. Enter information to connect to database and click on Connect Button




3. Choose one of three options to continue




4. Add New Visits




5. Look Up Patient




6. View/Edit Visits



7. Interview page: choose 1 of 2 options of THI and TFI interview



8. THI/TFI interview questions:



9. THI/TFI result page:

	Visit ID: 7	Date: Wed Nov	27 2019
Patient:	Adam Levine	THC #: 3	Visit no. 3
		estion 1: 4	
		estion 2: 4	
		estion 3: 4	
		estion 4: 4	
		estion 5: 4	
	() - V	estion 6: 4	
		estion 7: 4	
		estion 8: 4 estion 9: 4	
		estion 9: 4	
		estion 10: 4	
		stion 12: 4	
		stion 13: 4	
		stion 14: 4	
		stion 15: 4	
	Que	stion 16: 4	
	Que	stion 17: 2	
	Que	stion 18: 2	
		stion 19: 2	
		estion 20: 2	
		stion 21: 4	
		stion 22: 4	
		stion 23: 4	
	771-2-1-57	stion 24: 4	
		stion 25: 4	
	1. The state of th	nction score: 4	
		otion score: 3	
		strophic score	: 18
D		THI score: 92	. b dia
K	esult: GRADE 5	: Catastrophi	с папсісар
	CAN	CEL CAVE	
	CANO	CEL SAVE	

Project 3 - Team 7

TFI Result Back Visit ID: 7 Date: Wed Nov 27 2019 Patient: Adam Levine THC #: 3 Visit no. 3 Total INTRUSIVE score: 12 Total SENSE OF CONTROL score: 14 Total COGNITIVE score: 4 Total SLEEP score: 12 Total AUDITORY score: 13 Total RELAXATION score: 14 Total QUALITY OF LIFE score: 16 Total EMOTIONAL score: 14 Total TFI score: 45/100 Number of questions answered: 22 Result:

CANCEL

SAVE